


COLLOIDAL SILVER NANOPARTICLES

Product # SS-30



Colloidal Silver Nanoparticles
250 mL quantity shown



Transmission Electron Micrograph
37,500x magnification

Product Features

- Aqueous dispersion based on proprietary synthetic technology
- Highly monodisperse colloid, low polydispersity index (typically <0.5)
- High purity colloid for dependable and repeatable application results
- Available in quantities from 100 mL to 4 liters

Product Parameters

Material	Colloidal Silver Nanoparticles
Nominal Particle Size	28 nm
Average Zeta Potential	-40 mV
Form Factor	Colloid
Dispersion Solvent	Aqueous
Colloid Concentration	~15 ppm
Quantities Available	100 mL – 4 L

For pricing or additional product information, please contact our Sales Representative:

Phone: (540) 953-1785; E-mail: products@nanosonic.com

Products for research and development use only. Except where specifically enumerated by NanoSonic, Inc. by contractual agreement, all nanomaterials are for use in a research and development context only. They are not to be incorporated into commercialized food, household, cosmetic, or other applications where the buyer is not the end user. Further, materials sold by NanoSonic, Inc. are for the buying entity to use in an end-user context only. Only those individuals with proper technological background and training are permitted to work with NanoSonic nanomaterials, while using appropriate protective precautions. For further information, please consult the product MSDS, provided with each order and available upon request.

... Putting Nanotechnology To Work